LOGIKODIGIT

MA DM595 01 A_06-2022_ITA

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Venduto a _____

Via _____ **Località** _____

Joytech Healthcare Co., Ltd
No.365, Wuzhou Road, Yuhang Economic Development
Zone, Hangzhou City, 311100 Zhejiang P.R. China

Shanghai International Holding Corp. GmbH(Europe)
Eiffestrasse 80, 20537 Hamburg, Germany

Distribuito da: MORETTI S.P.A. Via Bruxelles, 3 - Melegnano 20022 Cavriaglia (Arezzo)
Tel. +39 055 96 21 11 - www.morettispa.com email: info@morettispa.com
MADE in P.R.C.

LOGIKODIGIT

MA DM595 01 A_06-2022_ENG

WIRST TYPE DIGITAL SPHYGMOMANOMETER

INSTRUCTION MANUAL

INDICE

1. CODE	PAG.4
2. INTRODUCTION.....	PAG.4
3. INTENDED USE	PAG.5
4. REGULATIONS AND DIRECTIVE OF REFERENCE	PAG.5
5. GENERAL WARNINGS	PAG.5
5.1 Battery Precautions	pag. 8
6. SYMBOLS.....	PAG.8
7. GENERAL DESCRIPTION.....	PAG.9
7.1 What is blood pressure?	pag. 9
7.2 Why measure your blood pressure?	pag. 9
8. CONTENTS	PAG.10
9. PARTS VIEW.....	PAG.11
9.1 Monitor unit	pag. 11
9.2 Display.....	pag. 11
9.2.1 Blood pressure WHO classification	pag. 12
10. IMPORTANT TESTING GUIDELINES	PAG.13
11. HOW TO USE.....	PAG.14
11.1. Quick start	pag. 14
11.2 Battery Installation	pag. 15
12. POWER ON.....	PAG.16
12.1 Power on	pag. 16
12.2 Positioning	pag. 16
12.3 Testing	pag. 17
12.4 Results.....	pag. 17
12.5 Irregular Heartbeat Indicator	pag. 17
12.6 Arm Shake Indicator.....	pag. 18
12.7 Power off.....	pag. 18
12.8 Average review of the last three measurements	pag. 18
12.9 Memory check.....	pag. 19

12.10 Memory deletion	pag. 19
12.11 Low Battery Indicator	pag. 19
12.12 Static Pressure Measurement.....	pag. 19
13. SYSTEM SETTINGS	PAG.20
13.1. Select memory group	pag. 20
13.2. Date/Time Settings	pag. 20
13.3. Time Format Settings	pag. 21
13.4. Unit settings	pag. 21
13.5. Volume Settings.....	pag. 22
13.6. Setting cuff positioning	pag. 22
13.7. Saving Settings	pag. 22
14. MAINTENANCE	PAG.22
15. CLEANING AND DISINFECTION	PAG.23
16. GENERAL CONDITIONS OF DISPOSAL	PAG.23
16.1 TREATMENT OF BATTERIES (Directive 2002/96/EC):.....	pag. 23
16.2 DISPOSAL OF WASTE BATTERIES (Directive 2006/66/EC):	pag. 24
17. TROUBLESHOOTING	PAG.24
18. Q&A.....	PAG.25
19. SPECIFICATIONS	PAG.26
20. ELECTROMAGNETIC COMPATIBILITY INFORMATION.....	PAG.28
21. WARRANTY	PAG.33
21.1 Repairing.....	pag. 33
21.2 Spare-parts	pag. 34
21.3 Exempt clauses	pag. 34

1. CODE

DM595 - wrist type Digital sphygmomanometer

2. INTRODUCTION

Thank you for purchasing LOGIKO sphygmomanometer by Moretti S.P.A. product range. This user's manual provides some suggestions as how to correctly use the product you have chosen and gives some valuable advice for your safety. Please read through the manual carefully before using the product. Should you have any queries, please contact your retailer for advice and assistance. In case of doubts, we recommend to contact the dealer who can help and advice you properly.

WARNING!

- Do not use the product for a purpose not indicated in this manual
- Moretti S.P.A declines all responsibilities for any consequences resulting from an incorrect use of this product and from unauthorized alteration to the frame of the product.
- The manufacturer reserves the right to change the information contained in this document without previous notice.
- This device is intended for non-invasive measuring of an adult individual's systolic, diastolic blood pressure and heart rate using the oscillometric method.
- Measurement position is on adult wrist only.

WARNING!

- The device is not intended to be used on infants and children.
- The device is designed for home or clinical use.
- All values can be read out in one LCD display.

3. INTENDED USE

The digital sphygmomanometer wrist type is a measurement device for systolic and diastolic blood pressure.

4. REGULATIONS AND DIRECTIVE OF REFERENCE

Safety standard:

EN 60601-1

Medical electrical equipment part 1: General requirements for safety

EMC standard:

EN 60601-1-2

Medical Electrical Equipment -- Part 1-2: General Requirements For Basic Safety And Essential Performance -- Collateral Standard: Electromagnetic Disturbances - Requirements And Tests

Performance standards:

IEC 80601-2-30

Medical electrical equipment - Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers.

EN 1060-3

Non-invasive sphygmomanometers - Supplementary requirements for electromechanical blood pressure measuring systems.

ISO 81060-2

Non-invasive sphygmomanometers - part 2: clinical validation of automated measurement type.

5. GENERAL WARNINGS

1. Do not confuse self-monitoring with self-diagnosis. Blood pressure measurements should only be interpreted by a health professional who is familiar with your medical history.
2. If you are taking medication, consult with your physician to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your physician.

3. For persons with irregular or unstable circulation resulting from diabetes, liver disease, arteriosclerosis or other medical conditions, there may be variations in blood pressure values measured at the wrist versus at the upper arm. Monitoring the trends in your blood pressure taken at either the arm or the wrist is nevertheless useful and important.
4. People suffering from vascular constriction, liver disorders or diabetes, people with cardiac pacemakers or a weak pulse, and women who are pregnant should consult their physician before measuring their blood pressure themselves. Different values may be obtained due to their condition.
5. People suffering from arrhythmias such as atrial or ventricular premature beats or atrial fibrillation only use this blood pressure monitor in consultation with your doctor. In certain cases oscillometric measurement method can produce incorrect readings.
6. Too frequent measurements can cause injury to the patient due to blood flow interference.
7. The cuff should not be applied over a wound as this can cause further injury.
8. DO NOT attach the cuff to a limb being used for IV infusions or any other intravascular access, therapy or an arterio-venous (A-V) shunt. The cuff inflation can temporarily block blood flow, potentially causing harm to the patient.
9. The cuff should not be placed on the arm on the side of a mastectomy. In the case of a double mastectomy use the side of the least dominant arm.
10. Pressurization of the cuff can temporarily cause loss of function of simultaneously used monitoring equipment on the same limb.
11. Check that operation of the unit does not result in prolonged impairment of the circulation of the patient.
12. Product is not intended for infants or individuals who cannot express their intentions.
13. Prolonged over-inflation of the bladder may cause ecchymoma of your arm.
14. Use only the approved arm cuff for this unit. Use of other arm cuffs may result in incorrect measurement results.

15. The system might produce incorrect readings if stored or used outside the manufacturer's specified temperature and humidity ranges.
16. Advising operator that Instruction manual/ Booklet must be consulted.
17. Do not use the device on vehicles that can influence measurement accuracy, such as patient transport in an ambulance or helicopter.
18. Contains small parts that may cause a choking hazard if swallowed by infants.

CAUTION!

Individuals with serious circulation problems may experience discomfort. Consult your physician prior to use.

Contact your physician if test results regularly indicate abnormal readings. Do not attempt to self-treat these symptoms without consulting your physician first.

Product is designed for its intended use only. Do not misuse in any way.

Product is not intended for infants or individuals who cannot express their intentions.

Do not disassemble or attempt to repair.

Do not use cell phones and other devices, which generate strong electrical or electromagnetic fields, near the device, as they may cause incorrect readings and interference or become interference source to the device.

Only use a recommended AC adaptor double-insulated complying with EN 60601-1 and EN 60601-1-2(see chapter 4). An unauthorized adaptor may cause fire and electric shock.

5.1 Battery Precautions

- Do not mix new and old batteries simultaneously.
- Replace batteries when Low Battery Indicator “” appears on screen.
- Be sure battery polarity is correct.
- Do not mix battery types. Long-life alkaline batteries are recommended.
- Remove batteries from device when not in operation for more than 3 months.
- Dispose batteries properly; observe local laws and regulations.

6. SYMBOLS

Caution

The product conforms to the requirements of the EC Directive MDD(93/42/EEC) on medical devices

Mandatory

Manufacturer

Prohibited

Authorised Representative in the European Community

Type BF Equipment

Keep off Sunlight

Instructions For Use
MUST be Consulted

Batch Lot

Serial Number

Manufacturing Date

Discard the used
product to the recycling
collection point according
to local regulations

IP22

Degree of protection
against dusts and liquids

7. GENERAL DESCRIPTION

7.1 What is blood pressure?

Pressure measures the action of the force exerted by blood on the walls of blood vessels. Systolic or maximum pressure indicates the maximum force exerted against the vessel walls during cardiac contraction. DIASTOLIC or minimum pressure indicates the force exerted on the vessels during cardiac relaxation. It is measured in mmHg (millimetres of mercury).


It is necessary to know that an individual's blood pressure fluctuates greatly and depends on many factors. It is generally lower in summer and higher in winter. Pressure can vary with atmospheric pressure and is greatly influenced by physical loads, emotional excitability, stress, meals, medication, alcohol, smoking, age, etc

It is recommended to note down your daily measurements and to consult with your physician who will tell you your normal blood pressure range.

NOTE: Blood pressure can change with the age, it is necessary to consult with a physician to know your normal blood pressure range. Never change drug dosage without consulting your physician.

7.2 Why measure your blood pressure?

Among the various health problems afflicting modern people, those associated with blood pressure are by far the most common. High blood pressure's dangerously strong correlation with cardiovascular disease and high morbidity has made measuring blood pressure as a necessity in identifying those at risk.

Example: Fluctuation within a day (Male 35 years old)

8. CONTENTS

1. Monitor unit

2. Owner's Manual


3. Plastic storage case


9. PARTS VIEW

9.1 Monitor unit


9.2 Display


9.2.1 Blood pressure WHO classification

The LOGIKO by Moretti monitors are equipped with a classification indicator based on established guidelines from the World Health Organization. The chart below (colour coded on monitor unit) indicates test results

■ Blood Pressure Classification Indicator

The following scheme has been developed by The World Health Organization (WHO), and it identifies the low- and high-risk blood pressure areas. This standard, however, is a general guideline as individuals' blood pressure varies among different people and different age groups, it is always subjective.

10. IMPORTANT TESTING GUIDELINES

1. Avoid eating, exercising, and bathing for 30 minutes prior to testing.
2. Sit in a calm environment for at least 5 minutes prior to testing.
3. Do not stand while testing. Sit in a relaxed position while keeping your arm level with your heart.
4. Avoid speaking or moving body parts while testing.
5. While testing, avoid strong electromagnetic interference such as microwave ovens and cell phones.
6. Wait 3 minutes or longer before re-testing.
7. Try to measure your blood pressure at the same time each day for consistency.
8. Test comparisons should only be made when monitor is used on the same arm, in the same position, and at the same time of day.
9. This blood pressure monitor is not recommended for people with severe arrhythmia.
10. Do not use this blood pressure monitor if the device is damaged.

11. HOW TO USE

11.1. Quick start

1. Install batteries (SEE FIGURE A)

Figure A

2. Remove any thick clothing from the wrist area (SEE FIGURE B)
3. Rest for a few minutes before the test. Wrap cuff around left wrist (See Figure C)

Figure B


Figure C

4. Sit in a comfortable position and place your wrist at heart level. see Figure D)
5. Press "START/STOP " to start the measurement. (See Figure E)

Figure D


Figure E

11.2 Battery Installation

Slide battery cover off as indicated by arrow.

Install 2 new AAA alkaline batteries according to polarity.


Close battery cover.

12. POWER ON

12.1 Power on


Press and hold "START/STOP" button to turn the unit on. The LCD screen will appear for one second as unit performs a quick diagnosis. A voice tone will indicate when unit is ready for testing.

Note: Unit will not function if residual air from previous testing is present in cuff. The LCD will flash "↓" until pressure is stabilized.

12.2 Positioning

Ensure that the wristband is level with the heart. When the wristband is not positioned correctly, the LCD screen displays POS and flashes as shown in the following interface:

By placing the wristband at the same level as the heart, the device will recognise the correct position.

12.3 Testing

After cuff inflation, air will slowly rise as indicated by the corresponding cuff pressure value. A flashing "♥" will appear simultaneously on screen signalling heart beat detection.

Note: Keep relaxed during testing. Avoid speaking or moving body parts.

12.4 Results

The screen will display systolic and diastolic blood pressure measurements. An indicator representing the current measurement will appear next to the corresponding WHO/WHO classification.

Note: Refer to the "Blood pressure WHO classification" section for detailed information.

12.5 Irregular Heartbeat Indicator

If the monitor detects an irregular heart rhythm two or more times during the measuring process, the Irregular Heartbeat Symbol "IHB" appears on screen along with measurement results. Irregular heartbeat rhythm is defined as rhythm that is either 25% slower or faster than the average rhythm detected while measuring systolic blood pressure and diastolic blood pressure. Consult

your physician if the Irregular Heartbeat Symbol "IHB" frequently appears with your test results.

12.6 Arm Shake Indicator

If there is arm movement during the measurement, the "⚡" icon may flash. Indicates that the measurement results may be inaccurate, and the situation will be recorded at the end of the measurement as a reminder.

12.7 Power off

The 'START/STOP' button can be pressed to switch off the unit in any mode. The unit will switch itself off after approximately 3 minutes of inactivity.

Safety precaution: If the pressure in the cuff becomes too high during the test, press the "START/STOP" button to switch off the unit. Cuff pressure dissipates quickly once the unit is switched off.

12.8 Average review of the last three measurements


With the power switched off, press the "MEM" button to activate the screen. After self-diagnosis, the screen will display the average results of the last 3 readings for the last group used. The "AVG" symbol will appear together with the corresponding WHO blood pressure indicator.

To check the average results of other groups, select the desired group before activating the "MEM" button

12.9 Memory check

By pressing the 'MEM' button again, previous results can be reviewed. After activating the review mode, you can press the 'MEM' button to scroll through all stored test results. The LCD will display the last memory as NO: 01.

Note: Previous results will only be displayed for the most recently used memory group. To check results stored in other memory groups, you must first select the desired group and then switch off the monitor.

12.10 Memory deletion

The memory of a selected group can be erased while in Review mode. Press and hold the 'START/STOP' button for about 3 seconds to delete all memory records of the selected group and then switch to test mode. Press and hold the 'START/STOP' button to turn off the unit.

Note: Memory cannot be recovered once deleted.

12.11 Low Battery Indicator


The unit will broadcast "Low Battery" when battery life is depleting and unable to inflate cuff for testing. The "⊗" appears simultaneously for approximately 5 seconds prior to shutting off. Replace batteries at this time. No memory loss will occur throughout this process.

12.12 Static Pressure Measurement

In the power down state, press and hold the "START/STOP" button, and then install the batteries. Until the LCD screen is full, release the "START/STOP" button.

When the LCD screen displays the double zero, the blood pressure metre is in static state.

Software version is displayed: 10 is a software version in the figure.

Note: Only Service personnel can have access to this mode.

13. SYSTEM SETTINGS

With the power off, press and hold the 'START/STOP' button for 3 seconds to activate the system settings. The memory group icon flashes.

13.1. Select memory group

In the system setting mode, measurement data can be stored in 2 different groups. This allows multiple users to save individual test results (up to 150 memories per group). Press the 'MEM' button to choose a group setting. Test results will automatically be stored in each selected group.

13.2. Date/Time Settings

Press the 'START/STOP' button again to set the time/date mode. Set the year first by adjusting the 'MEM' button. Press the "START/STOP" button again to confirm the current month. Continue setting the day, hour and minute in the same way. Each time the 'START/STOP' button is pressed, the selection is stopped and continued in succession (month, day, hour, minute, 12/24 hour)

13.3. Time Format Settings

Press the 'START/STOP' button again to set the time format mode. Set the time format by adjusting the "MEM" button. EU means European time, US means American time.

13.4. Unit settings

Press the 'START/STOP' button again to set the unit. Press the "MEM" button to change settings.

13.5. Volume Settings

Press the 'START/STOP' button to enter the volume setting mode. Set the voice volume by adjusting the 'MEM' button. Smaller " " is for lower volume. There are six volume levels.

13.6. Setting cuff positioning

Press the 'START/STOP' button to enter setting. Set the cuff positioning function at heart level to ON or OFF by pressing the "MEM" button.

13.7. Saving Settings

While in any setting mode, press and hold the "START/STOP" button to switch off the unit. All information will be saved.

Note: The unit will automatically save all information and shut down if left idle for 3 minutes.

14. MAINTENANCE

The LOGIKO devices by Moretti S.P.A are checked carefully and supplied with EC mark, once launched on the market.

For patient and physician's safety we recommend to let check the device by the manufacturer or by laboratory authorized every 2 years. In case of repair, use only original spare parts and accessories.

15. CLEANING AND DISINFECTION

- DO NOT immerse the device!
- Cleaning: use only a damp cloth. Pay attention to the display
- Never use acids, alkalis or solvents such as acetone
- Avoid dropping, slamming, or throwing the unit
- Avoid extreme temperatures. Do not use outdoors
- Do not disassemble product
- It is recommended the performance should be checked every 2 years
- Remove batteries when not in operation for an extended period of time

16. GENERAL CONDITIONS OF DISPOSAL

In case of disposal don't use the inserting container for municipal waste. We recommend to dispose the sphygmomanometer in the appropriate disposal areas for recycling.

16.1 TREATMENT OF BATTERIES (Directive 2002/96/EC):

At the end of its life, the product must not be disposed of along with other domestic waste. The users must dispose of this equipment by bringing it to a specific recycling point for electric and electronic equipment or at retailers that provide this service. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. To remark the need to dispose of electrical equipment separately, the products is marked with crossed mobile waste bin.

**16.2 DISPOSAL OF WASTE BATTERIES (Directive 2006/66/EC):**

This symbol on the battery or on the packaging indicates that the battery provided with this product shall not be treated as household waste. By ensuring these batteries are disposed of correctly, you will help prevent potentially negative consequences for the environment and human health which could otherwise be caused by inappropriate waste handling of the battery. The recycling of the materials will help to conserve natural resources. At the end of their life, hand the batteries over to the applicable collection points for the recycling of waste batteries. For more detailed information about recycling of this product or battery, please contact your local Civic Office, your household waste disposal service or the shop where you purchased the product.

17. TROUBLESHOOTING

Abnormal phenomenon	Cause analysis	Processing method
Abnormal sphygmomanometer	The armband is tied too tight or too loose, Or the arm strap is tied incorrectly;	Roll the armband correctly
	Move the arm during measurement or Electronic sphygmomanometer	Stay quiet, keep your arm steady, and do not move the monitor
	Speaking, nervous or emotional during measurement	Instead of talking, take deep breaths to calm your mood and relax your body
	Incorrect measurement posture	Adjust posture, see "Blood pressure gauge Wearing"
	There is interference in charging process or improper operation in measuring process	See operation Instructions.

The following table shows the error signs that may occur during measurement, possible causes and handling methods. Please measure again using the correct method

Error	The cause of the problem	The solution
Er1	Can't detect high and low pressure	Please fasten the cuff before measuring
Er2	Cuff too loose or loose	Please fasten the cuff before measuring
Er3	Improper compression caused by arm or body movement	Hold the arm or body still and measure again
Er4	The pressure exceeds 300mmHg	Please fasten the cuff before measuring
Er5	The pressure exceeds 15mmHg for 3 minutes	Check whether the cuff is knotted or the vent valve is blocked. If the
	Battery dead	Replace the battery or connect the power adapter (if any).

Note: If you cannot solve the abnormal situation by yourself, you can consult the manufacturer or the manufacturer's designated unit by phone. It is forbidden to disassemble and repair without permission. If necessary, professional maintenance personnel can ask the manufacturer for the list of components and circuit schematic diagram.

18. Q&A

- Q:** What is the difference between measuring blood pressure at home or at a professional healthcare clinic?
- A:** Blood pressure readings taken at home are now seen to give a more accurate account as they better reflect your daily life. Readings can be elevated when taken in a clinical or medical environment. This is known as White Coat Hypertension and may be caused by feeling anxious or nervous.

Note: Abnormal test results may be caused by:

1. Improper cuff placement
Make sure cuff is snug-not too tight or too loose.
2. Improper body position
Make sure to keep your body in an upright position.
3. Feeling anxious or nervous
Take 2-3 deep breaths, wait a few minutes and resume testing.

Q: What causes different readings?

A: Blood pressure varies throughout the course of a day. Many factors may affect an individual's blood pressure.

Q: Should I apply the cuff to the left or right arm? What is the difference?

A: Either arm can be used when testing, however, when comparing results, the same arm should be used. Testing on your left arm may provide more accurate results as it is located closer to your heart.

Q: What is the best time of day for testing?

A: Morning time or any time you feel relaxed and stress free.

19. SPECIFICATIONS

Description	Wrist-type Fully Automatic Blood Pressure Monitor
Model	DM595
Display	LCD Digital Display Size:32mm x 40.5mm(1.26" x 1.59")
Measurement Method	Oscillometric Method

Pressure measurement range	Systolic	60mmHg~260 mmHg
	Diastolic	30mmHg~200 mmHg
	Pressure	0mmHg~299mmHg
	Accuracy	± 3mmHg
	Pulse	30 ~ 180 Beats/Minute
	Accuracy	± 5%
Pressurization	Automatic Pressurization	
Memory	2x150 Memories in Tow Groups with Date and Time	
Functions	Irregular Heartbeat Detection	
	WHO Classification Indicator	
	Last 3 Results Average	
	Low Battery Detection	
	Automatic Power-Off	
	Backlight	
Power Source	2 Alkaline Batteries Size AAA	
Battery Life	Approximately 2 months at 3 tests per day	
Unit Weight	Approx. 67g (2.36oz) (Excluding Battery)	
Unit Dimensions	Approx. 84mm x 62mm x 23.6mm(L x W x H) (3.31" x 2.44" x 0.93")	
Cuff Circumference	Fits wrist circumference 13.5-21.5 cm(5.3"-8.5")	
Operating Environment	Temperature	10°C ~ 40°C (50°F~104°F)
	Humidity	15%~93%RH
	Pressure	800hPa~1060hPa
Storage Environment	Temperature	-25°C~55°C (-13°F~131°F)
	Humidity	≤93% RH
Ingress Protection Rating	IP 22	
Classification	Internal Powered Equipment Type BF  Cuff is the Applied Part	

Specifications are subject to change without notice.

20. ELECTROMAGNETIC COMPATIBILITY INFORMATION

The device satisfies the EMC requirements of the international standard IEC 60601-1-2. The requirements are satisfied under the conditions described in the table below. The device is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect

the device. Use of the unit in conjunction with non-approved accessories can affect the device negatively and alter the electromagnetic compatibility. The device should not be used directly adjacent to or between other electrical equipment.

Table 1

Guidance and declaration of manufacturer-electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment -guidance
Radiated emission CISPR 11	Group 1, ClassB	The device uses RF energy only for its internal function. Therefore, its emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted emission CISPR 11	N/A	
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	N/A	

Table 2

Guidance and declaration of manufacturer-electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient/burst IEC 61000-4-4	± 2 kV , 100kHz, for AC power port	N/A	
Surge IEC 61000-4-5	± 0.5kV, ± 1kV (differenziale)	N/A	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0° 0 % UT; 250/300 cycle	N/A	
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m; 50Hz or 60Hz	30 A/m; 50Hz or 60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Radiated RF EM fields IEC 61000-4-3	3V/m or 10 V/m 80MHz-2.7 Ghz 80%AM at 1kHz	3V/m or 10 V/m 80MHz-2.7 Ghz 80%AM at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance 80 MHz to 800 MHz 800 MHz to 2.7 Ghz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted disturbances Induced by RF fields IEC 61000-4-6	3 V in 0.15 MHz-80 MHz 6 V in ISM and/or amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz	3 V in 0.15 MHz-80 MHz 6 V in ISM and/or amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance 80 MHz to 800 MHz 800 MHz to 2.7 Ghz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

Table 3

[illegible]

2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n,RFID 2450,LTE Band 7	Pulse Modulation 217Hz	2	0.3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse Modulation 217Hz	0.2	0.3	9
5500						
5785						

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device		
The device is intended for use in an electromagnetic environment in which radiated therefore disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.		
Rated maximum output power of transmitter	Separation distance according to frequency of transmitter	
	80 MHz to 800 MHz $d=[3.5/E_1]\sqrt{P}$	80 MHz to 2.7 GHz $d=[7/E_1]\sqrt{P}$
0.01	0.12	0.23
0.1	0.38	0.73
1	1.2	2.3
10	3.8	7.3
100	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.		
NOTE1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.		
NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

21. WARRANTY

Moretti products are guaranteed from material or manufacturing faults for 2 years from the purchasing date, except possible exclusion or restriction as follows. The warranty shall not be applied in the possible damages caused by improper use, abuse or alteration, and the warranty will not be valid if the instructions for use are not strictly adhered to. The correct intended use is specified in this manual. Moretti is not responsible for consequent damages, personal injuries or whatever caused by or in relation to wrong installation or improper use. Moretti warranty does not cover damages resulting from: natural disaster, not authorized maintenance or repairs, faults caused by problems on electricity supply (when necessary), use of spare parts not covered by Moretti, improper use, not authorized alteration, shipment damages (different from original Moretti shipment), or in case of insufficient maintenance as indicated in the manual. The warranty doesn't cover components subject to wear and tear during the correct use of the device.

21.1 Repairing

-Warranty repair-

If a Moretti item presents material or manufacturing faults during the warranty period, Moretti will confirm with customer if the fault can be covered from warranty. Moretti, at its unquestionable discretion, can repair or replace the item, by a Moretti dealer or to Moretti headquarters. Labor cost can be charged to Moretti if the repair is covered by warranty. A repair or a replacement doesn't extend the warranty.

-Repair not covered by warranty-

A product out of warranty can be sent after Moretti authorization. The labor and shipping costs for good out of warranty are to be paid by the customer or by the dealer. The repairs are guaranteed for 6 months from the good received.

**-Non-defective products-**

The customer will be notified if, after examining and tried a returned product, Moretti concludes that the product is not defective. The product will be returned to the customer and will be liable to pay the costs caused by the return.

21.2 Spare-parts

Moretti's genuine parts are guaranteed for 6 months with effect from the day you receive the replacement.

21.3 Exempt clauses

Moretti does not offer any other declaration, explicit or implicit warranty or conditions, including possible declarations, warranties or conditions of merchantability, fitness for a specific purpose, non infringement and non interference, all but what expressly specified in this warranty. Moretti does not guarantee the nonstop and faultless usage. The duration of possible implicit warranties which can be imposed by the law is limited by the warranty period, in the limits of law. Some states or countries don't allow limiting the implicit warranty or the exclusion or the limitation for accidental damages. In such countries, some of those exclusions or limitation may not be applied to the user. The present warranty may be modified without prior notification.

LOGIKODIGIT

WARRANTY CERTIFICATE

Product _____

Purchasing date _____

Authorized dealer _____

Street _____ Place _____

Sold to _____

Street _____ Place _____

Joytech Healthcare Co., Ltd
No.365, Wuzhou Road, Yuhang Economic Development
Zone, Hangzhou City, 311100 Zhejiang P.R. China

Shanghai International Holding Corp. GmbH(Europe)
Eiffestrasse 80, 20537 Hamburg, Germany

Distributed by: MORETTI S.P.A. Via Bruxelles, 3 - Melegnano 20022 Cavriaglia (Arezzo)
Tel. +39 055 96 21 11 - www.morettispa.com email: info@morettispa.com
MADE in P.R.C.

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MORETTI S.p.A.