GIMA

Test Monofase Cocaina Strisce Urina
One Step Cocaine Test Strip (Urine)
1 Etape Test de Cocaine sur Bandelette (Urine)
Einstufen Kokain Teststreifen (Urin)
Prueba de Cocaína en Un Solo Paso en Tira (Orina)
Tira para Teste de Cocaína Em Um Só Passo (Urina)
Τεστ Κokaίνης Ταινίας Ενός Σταδίου (Ούρων)
فحص أحادي الطور كوكاينا شرائح بول

MANUALE D'USO
OPERATOR'S MANUAL
MANUEL D'UTILIZATION
BEDIENUNGSANLEITUNG
MANUAL DE USO
MANUAL DE USO
Εγχειρίδιο χρήσης
دليل للإرشادات

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
BIBLIOGRAFIA

1. Stewart DI, Inoba T, Ducassen M, Kalow W. Clin. Pharmacol. Ther. 1979; 25: 264
2. Ambre J. J. Anal. Toxicol. 1985; 9: 241
3. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982: 488
4. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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 Fabbricante GIMA Spa
 Via Monza, 102 - 20060
 Gessate (MI) - Italia

One Step Cocaine Test Strip (Urine)

A rapid, one step test for the qualitative detection of Cocaine metabolites in human urine. For medical and other professional in vitro diagnostic use only.

INTENDED USE

The COC One Step Cocaine Test Strip (Urine) is a rapid chromatographic immunoassay for the qualitative detection of Cocaine metabolite, Benzoylcegonine, in human urine at a cut-off concentration of 300 ng/mL. This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert. **This assay provides only a qualitative, preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are used.**

SUMMARY

Cocaine, is a potent central nervous system (CNS) stimulant and a local anesthetic. Initially, it brings about extreme energy and restlessness while gradually resulting in tremors, over-sensitivity and spasms. In large amounts, Cocaine causes fever, unresponsiveness, and difficulty in breathing and unconsciousness. Cocaine is often self-administered by nasal inhalation, intravenous injection and free-base smoking. It is excreted in the urine in a short time primarily as Benzoylcegonine.^{1,2} Benzoylcegonine, a major metabolite of Cocaine, has a longer biological half-life (5 -8 hours) than Cocaine (0.5 -1.5 hours), and can generally be detected for 24-48 hours after Cocaine exposure.² The COC One Step Cocaine Test Strip (Urine) is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes a monoclonal antibody to selectively detect elevated levels of Cocaine metabolite in urine. The COC One Step Cocaine Test Strip (Urine) yields a positive result when the Cocaine metabolite in urine exceeds 300 ng/mL. This is the suggested screening cut-off for positive specimens set by the Substance Abuse and Mental Health Services Administration (SAMHSA, USA).

PRINCIPLE

The COC One Step Cocaine Test Strip (Urine) is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against the drug conjugate for binding sites on the antibody. During testing, a urine specimen migrates upward by capillary action. Benzoylcegonine, if present in the urine specimen below 300 ng/mL, will not saturate the binding sites of antibody in the test. The antibody coated particles will then be captured by immobilized Benzoylcegonine conjugate and a visible colored line will appear in the test line region. The colored line will not form in the test line region if the Benzoylcegonine level is above 300 ng/mL because it will saturate all the binding sites of antibodies. A drug-positive urine specimen will not generate a colored line in the test line region because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The test contains mouse monoclonal anti-Benzoylcegonine antibody-coupled particles and Benzoylcegonine-protein conjugate. A goat antibody is employed in the control line system.

PRECAUTIONS

- For medical and other professional in vitro diagnostic use only. Do not use after the expiration date.
- The test should remain in the sealed pouch until use.
- All specimens must be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test should be discarded according to local regulations.

STORAGE AND STABILITY

Store as packaged in the sealed pouch either at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. **DO NOT FREEZE**. Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.

Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to testing. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.

MATERIALS

Materials Provided

- Test strips
- Package insert


Materials Required But Not Provided


- Specimen collection container
- Timer

DIRECTIONS FOR USE

Allow the test, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

1. Bring the pouch to room temperature before opening it. Remove the test strip from the sealed pouch and use it as soon as possible.
2. With arrows pointing toward the urine specimen, immerse the test strip vertically in the urine specimen for at least 10-15 seconds. Do not pass the maximum line (MAX) on the test strip when immersing it. See the illustration below.
3. Place the test strip on a non-absorbent flat surface, start the timer and wait for the colored line(s) to appear. Read results at 5 minutes. Do not interpret the result after 10 minutes.

INTERPRETATION OF RESULTS

(Please refer to the illustration above)

NEGATIVE:* Two lines appear. One colored line should be in the control line region (C), and another apparent colored line should be in the test line region (T). This negative result indicates that the Benzoylcegonine concentration is below the detectable level (300 ng/mL).

*NOTE: The shade of color in the test line region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

POSITIVE: One colored line appears in the control line region (C). No line appears in the test line region. This positive result indicates that the Benzoylcegonine concentration is above the detectable level (300 ng/mL).

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A colored line appearing in the control line region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATION

1. The COC One Step Cocaine Test Strip (Urine) provides only a qualitative, preliminary analytical result. A secondary quantitative analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.^{3,4}
2. It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
3. Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
4. A positive result does not indicate level of intoxication, administration route or concentration in urine.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
6. Test does not distinguish between drugs of abuse and certain medications.

PERFORMANCE CHARACTERISTICS

Accuracy

A side-by-side comparison was conducted by laboratory personnel using the COC One Step Cocaine Test Strip (Urine) and a commercially available rapid test. Testing was performed on 300 clinical specimens previously collected from subjects present for Drug Screen Testing. Ten percent of the specimens employed were either at -25% or +25% level of the cut-off concentration of 300 ng/mL Benzoylcegonine. Presumptive positive results were confirmed by GC/MS. The following results were tabulated:

Method	Other COC Rapid Test		Total Results
	Results	Positive	Negative
COC One Step Test Strip	Positive	136	0
	Negative	7	157
Total Results		143	157
% Agreement		95%	>99%
			98%

When compared to GC/MS at the cut-off of 300 ng/mL, the following results were tabulated:

COC One Step Test Strip	Method	GC/MS		Total Results
	Results	Positive	Negative	
	Positive	119	17	136
	Negative	5	159	164
	Total Results	124	176	300
	% Agreement	96%	90%	93%

Analytical Sensitivity

A drug-free urine pool was spiked with Benzoylcegonine at the following concentrations: 0 ng/mL, 150 ng/mL, 225 ng/mL, 300 ng/mL, 375 ng/mL and 450 ng/mL. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

Benzoylcegonine Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
150	-50%	30	30	0
225	-25%	30	30	0
300	Cut-off	30	4	26
375	+25%	30	0	30
450	+50%	30	0	30

Analytical Specificity

The following table lists compounds that are positively detected in urine by the COC One Step Cocaine Test Strip (Urine) at 5 minutes.

Compound	Concentration (ng/mL)
Benzoylcegonine	300
Cocaine HCl	780
Cocaethylene	12,500
Ecgonine HCl	32,000

Precision

A study was conducted at three physicians' offices by untrained operators using three different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to GC/MS, no Benzoylcegonine, 25% Benzoylcegonine above and below the cut-off and 50% Benzoylcegonine above and below the 300 ng/mL cut-off was provided to each site. The following results were tabulated:

Benzoylcegonine Concentration (ng/mL)	n per Site	Site A		Site B		Site C	
		Pos	Neg	Pos	Neg	Pos	Neg
0	15	0	14 ^a	0	15	0	15
150	15	1	14	0	15	1	14
225	15	11	4	10	5	7	8
375	15	15	0	15	0	15	0
450	15	15	0	15	0	14	1
Non Valid	15	16/16		15/15		15/15	

*Note: Non-valid results were obtained in this treatment. Non-valid tests were provided as part of this study to ensure that readers would accurately identify non-valid test results.

Effect of Urinary Specific Gravity

Fifteen urine specimens of normal, high, and low specific gravity ranges were spiked with 150 ng/mL and 450 ng/mL of Benzoylcegonine. The COC One Step Cocaine Test Strip (Urine) was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Benzoylcegonine to 150 ng/mL and 450 ng/mL. The spiked, pH-adjusted urine was tested with the COC One Step Cocaine Test Strip (Urine) in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-negative urine or Benzoylcegonine positive urine. The following compounds show no interference when tested with the COC One Step Cocaine Test Strip (Urine) at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds

Acetaminophen	Diazepam	Methadone	Prednisone
Acetophenetidin	Diclofenac	Methoxyphenamine	Procaine
N-Acetylprocainamide	Diflunisal	(±)3,4-Methylenedioxy- amphetamine	Promazine
Acetylsalicylic acid	Digoxin	(±)-3,4-Methylenedioxy- methamphetamine	Promethazine
Aminopyrine	Diphenhydramine	Morphine-3-β -D glucuronide	D,L-Propranolol
Amitriptyline	Doxylamine	Morphine Sulfate	D-Propoxyphene
Amobarbital	Ecgonine methylester	Naloxone	D-Pseudoephedrine
Amoxicillin	(-)-℥ -Ephedrine	Naltrexone	Quinidine
Ampicillin	Erythromycin	Naproxen	Quinine
L-Ascorbic acid	β-Estradiol	Niacinamide	Ranitidine
D,L-Amphetamine sulfate	Estrone-3-sulfate	Nifedipine	Salicylic acid
Apomorphine	Ethyl-p-aminobenzoate	Norcodein	Secobarbital
Aspartame	Fenopropfen	Norethindrone	Serotonin
Atropine	Furosemide	D-Norpropoxyphene	Sulfamethazine
Benzilic acid	Gentisic acid	Noscapine	Sulindac
Benzoic acid	Hemoglobin	Oxalic acid	Temazepam
Benzphetamine	Hydralazine	Oxazepam	Tetracycline
Bilirubin	Hydrochlorothiazide	Oxolinic acid	Tetrahydrocortisone, 3-acetate
(±)-Brompheniramine	Hydrocodone	Oxycodone	Tetrahydrocortisone 3-(β-D glucuronide)
Caffeine	Hydrocortisone	Papaverine	Tetrahydrozoline
Cannabidiol	O-Hydroxyhippuric acid	Penicillin-G	Thebaine
Cannabinol	p-Hydroxy- methamphetamine	Pentobarbital	Thiamine
Chloralhydrate	3-Hydroxytyramine	Perphenazine	Thioridazine
Chloramphenicol	Ibuprofen	Phencyclidine	D,L-Tyrosine
Chlordiazepoxide	Imipramine	Phenelzine	Tolbutamide
Chlorothiazide	Iproniazid	Phenobarbital	Triamterene
(±)-Chlorpheniramine	(±)-Isoproterenol	Phentermine	Trifluoperazine
Chlorpromazine	Isoxsuprine	L-Phenylephrine	Trimethoprim
Chlorquine	Ketamine	β-Phenylethylamine	Trimipramine
Cholesterol	Ketoprofen	Phenylpropanolamine	Tryptamine
Clomipramine	Labetalol	Prednisolone	D,L-Tryptophan
Clonidine	Levorphanol		Tyramine
Codeine	Loperamide		Uric acid
Cortisone	Maprotiline		Verapamil
(-)Cotinine	Meperidine		Zomepirac
Creatinine	Meprobamate		
Deoxycorticosterone			
Dextromethorphan			

BIBLIOGRAPHY

1. Stewart DI, T Inoba, M Ducassen, W Kalow. Clin. Pharmacol. Ther. 1979; 25: 264
2. Ambre J.J. Anal. Toxicol. 1985; 9: 241
3. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982: 488
4. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

Index of Symbols

	Attention, see instruction for use		Test per kit		Manufacturer
	For <i>in vitro</i> diagnostic use only		Use by		Do not reuse
	Store between 2-30° C		Lot Number		Catalog 24540

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Cholesterol
Clomipramine
Clonidine
Codeine
Cortisone
(-)Cotinine
Creatinine
Deoxycorticosterone
Dextromethorphan

Ketamine
Ketoprofen
Labetalol
Levorphanol
Loperamide
Maprotiline
Meperidine
Meprobamate

Phencyclidine
Phenelzine
Phenobarbital
Phentermine
L-Phenylephrine
β-Phenylethylamine
Phenylpropanolamine
Prednisolone

Trimethoprim
Trimipramine
Tryptamine
D,L-Tryptophan
Tyramine
Uric acid
Verapamil
Zomepirac

BIBLIOGRAPHIE

1. Stewart DJ, T Inoba, M Ducassen, W Kalow. Clin. Pharmacol. Ther. 1979; 25: 264
2. Ambre J. J. Anal. Toxicol. 1985; 9: 241
3. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982: 488
4. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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Acetaminophen	Diazepam	Methadone	Prednisone
Acetophenetidin	Diclofenac	Methoxyphenamine	Procaine
N-Acetylprocainamide	Diflunisal	(±)3,4-Methylenedioxy- amphetamine	Promazine
Acetylsalicylic acid	Digoxin	(±)-3,4-Methylenedioxy- methamphetamine	Promethazine
Aminopyrine	Diphenhydramine	Morphine-3-β -D glucuronide	D,L-Propranolol
Amitriptyline	Doxylamine	Morphine Sulfate	D-Propoxyphene
Amobarbital	Ecgonine methylester	Nalidixic acid	D-Pseudoephedrine
Amoxicillin	(-)-N-Ephedrine	Naloxone	Quinidine
Ampicillin	Erythromycin	Naltrexone	Quinine
L-Ascorbic acid	β-Estradiol	Naproxen	Ranitidine
D,L-Amphetamine sulfate	Estrone-3-sulfate	Niacinamide	Salicylic acid
Apomorphine	Ethyl-p-aminobenzoate	Nifedipine	Secobarbital
Aspartame	Fenopropfen	Norcodein	Serotonin
Atropine	Furosemide	Norethindrone	Sulfamethazine
Benzilic acid	Gentisic acid	D-Norpropoxyphene	Sulindac
Benzoic acid	Hemoglobin	Noscapine	Temazepam
Benzphetamine	Hydralazine	Oxalic acid	Tetracycline
Bilirubin	Hydrochlorothiazide	Oxazepam	Tetrahydrocortisone, 3-acetate
(±)-Brompheniramine	Hydrocodone	Oxolinic acid	Tetrahydrocortisone 3-(β-D glucuronide)
Caffeine	Hydrocortisone	Oxycodone	Tetrahydrozoline
Cannabidiol	O-Hydroxyhippuric acid	Papaverine	Thebaine
Cannabinol	p-Hydroxy-	Penicillin-G	Thiamine
Chloralhydrate	Methamphetamine	Pentobarbital	Thioridazine
Chloramphenicol	3-Hydroxytyramine	Perphenazine	D,L-Tyrosine
Chlordiazepoxide	Ibuprofen	Phencyclidine	Tolbutamide
Chlorothiazide	Imipramine	Phenelzine	Triamterene
(±)-Chlorpheniramine	Iproniazid	Phenobarbital	Trifluoperazine
Chlorpromazine	(±)-Isoproterenol	Phentermine	Trimethoprim
Chlorquine	Isoxsuprine	L-Phenylephrine	Trimipramine
Cholesterol	Ketamine	β-Phenylethylamine	Tryptamine
Clomipramine	Ketoprofen	Phenylpropanolamine	D,L-Tryptophan
Clonidine	Labetalol	Prednisolone	Tyramine
Codeine	Levorphanol		Uric acid
Cortisone	Loperamide		Verapamil
(-)-Cotinine	Maprotiline		Zomepirac
Creatinine	Meperidine		
Deoxycorticosterone	Meprobamate		
Dextromethorphan			

BIBLIOGRAFIA

1. Stewart DI, T Inoba, M Ducassen, W Kalow. Clin. Pharmacol. Ther. 1979; 25: 264
2. Ambre J. J. Anal. Toxicol. 1985; 9: 241
3. Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982: 488
4. Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986

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Non Cross-Reacting Compounds

Acetaminophen	Diazepam	Methadone	Prednisone
Acetophenetidin	Diclofenac	Methoxyphenamine	Procaine
N-Acetylprocainamide	Diflunisal	(±)3,4-Methylenedioxy-amphetamine	Promazine
Acetylsalicylic acid	Digoxin	(±)-3,4-Methylenedioxy-methamphetamine	Promethazine
Aminopyrine	Diphenhydramine	Morphine-3-β-D-glucuronide	D,L-Propranolol
Amitriptyline	Doxylamine	Morphine Sulfate	D-Propoxyphene
Amobarbital	Ecgonine methylester	Nalidixic acid	D-Pseudoephedrine
Amoxicillin	(-)-Ψ-Ephedrine	Naloxone	Quinidine
Ampicillin	Erythromycin	Naltrexone	Quinine
L-Ascorbic acid	β-Estradiol	Naproxen	Ranitidine
D,L-Amphetamine sulfate	Estrone-3-sulfate	Niacinamide	Salicylic acid
Apomorphine	Ethyl-p-aminobenzoate	Nifedipine	Secobarbital
Aspartame	Fenopropfen	Norcodein	Serotonin
Atropine	Furosemide	Norethindrone	Sulfamethazine
Benzilic acid	Gentisic acid	D-Norpropoxyphene	Sulindac
Benzoic acid	Hemoglobin	Noscapine	Temazepam
Benzphetamine	Hydralazine	Oxalic acid	Tetracycline
Bilirubin	Hydrochlorothiazide	Oxazepam	Tetrahydrocortisone, 3-acetate
(±)-Brompheniramine	Hydrocodone	Oxolinic acid	Tetrahydrocortisone 3-(β-D glucuronide)
Caffeine	Hydrocortisone	Oxycodone	Tetrahydrozoline
Cannabidiol	O-Hydroxyhippuric acid	Papaverine	Thebaine
Cannabinol	p-Hydroxy-	Penicillin-G	Thiamine
Chloralhydrate	Methamphetamine	Pentobarbital	Thioridazine
Chloramphenicol	3-Hydroxytyramine	Perphenazine	D,L-Tyrosine
Chlordiazepoxide	Ibuprofen	Phencyclidine	Tolbutamide
Chlorothiazide	Imipramine	Phenelzine	Triamterene
(±)-Chlorpheniramine	Iproniazid	Phenobarbital	Trifluoperazine
Chlorpromazine	(±)-Isoproterenol	Phentermine	Trimethoprim
Chlorquine	Isoxsuprine	L-Phenylephrine	Trimipramine
Cholesterol	Ketamine	β-Phenylethylamine	Tryptamine
Clomipramine	Ketoprofen	Phenylpropanolamine	D,L-Tryptophan
Clonidine	Labetalol	Prednisolone	Tyramine
Codeine	Levorphanol		Uric acid
Cortisone	Loperamide		Verapamil
(-)-Cotinine	Maprotiline		Zomepirac
Creatinine	Meperidine		
Deoxycorticosterone	Meprobamate		
Dextromethorphan			

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