GIMA

Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
Italia: tel. 199 400 401 - fax 199 400 403 Export: tel. +39 02
953854209/221/225 fax +39 02 95380056 gima@gimaitaly.
com - export@gimaitaly.com www.gimaitaly.com

PROFESSIONAL MEDICAL PRODUCTS

**TEST PER LA VISIONE STEREOSCOPICA – TEST FOR
STEREOSCOPIC VISION – EXAMEN DE LA VISION
STÉRÉOSCOPIQUE – TEST FÜR STEREOSKOPISCHES SEHEN –
TEST PARA LA VISIÓN ESTEREOSCÓPICA – TESTE PARA A VISÃO
ESTEREOSCÓPICA – ΤΕΣΤ ΣΤΕΡΕΟΣΚΟΠΙΚΗΣ ΟΡΑΣΗΣ –**

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- Rowe FJ: *Clinical orthoptics*. Wiley-Blackwell 3^a ediz. 2012

- Duckman RH: *Visual development, diagnosis, and treatment of the pediatric patient*. Lippincott Williams & Wilkins 2006

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ENGLISH

Evaluates the cooperation between both eyes in forming a single image, starting from the two images that are supplied separately to the brain by the two eyes.

Can be used for children over one year of age and adults.

METHOD AND HOW TO USE:

The test consists of a card showing 4 shapes:

a ball, a fish, a star and a tree.

Familiar images have been used, which can be recognised and named both by adults and by very young children.

The person being tested is asked to indicate what he/she sees in the test and his/her eye movements are observed.

The ball is the only element visible using just one eye. The other shapes (each of which is placed at different depth levels) are perceived only by someone with binocular vision and visual acuity over 2-3/10. For adults and cooperative children the card should be held at a distance of 30-40 cm, with a request to indicate the shapes seen in it. If only the ball is seen, the result is pathological. If all four shapes can be identified, then the result is normal.

Briefly:

Positive result:

All objects are identified and named

Unclear result:

Only one of the images is recognised

Negative result:

None of the three-dimensional images is recognised.

INDICATIONS FOR EVALUATION OF THE TEST:

Suspected strabismus, amblyopia or reduced sight in one or both eyes or as a screening test for stereoscopic vision.

If a child does not cooperate, the test can be shown and the tester can observe whether the child tries to grasp one of the objects apart from the ball (the fish, the star or the tree). If this occurs, the test is normal.

If it does not happen, the test cannot be considered pathological for certain and may be repeated some time later.

Important: in the event of a negative or unclear result, it is advisable to carry out a complete eye examination.

GIMA WARRANTY CONDITIONS

Congratulations for purchasing a GIMA product. This product meets high qualitative standards both as regards the material and the production. The warranty is valid for 12 months from the date of supply of GIMA. During the period of validity of the warranty, GIMA will repair and/or replace free of charge all the defected parts due to production reasons. Labor costs and personnel traveling expenses and packaging not included. All components subject to wear are not included in the