medel[®]

CARDI WATCH

MANUALE D'USO
USER MANUAL

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+39 02 83451193

1. PRODUCT SPECIFICATIONS

Thank you for purchasing the Medel CARDIO WATCH.

Check the external integrity of the package and the contents. Before use, make sure that the appliance and accessories do not show any obvious damage and that the packaging material has been removed. If in doubt, do not use the appliance and consult your dealer or contact the customer Service indicated.

PACKAGE CONTENTS:

- 1 bracelet-shaped activity sensor
- 1 x USB charging cable

Why use an activity sensor?

Movement is the fundamental condition for a healthy life. Physical activity helps keep fit and reduces the risk of developing cardiovascular system diseases. Physical activity means any movement that involves energy consumption. Daily activities such as walking, cycling, climbing stairs and doing housework contribute greatly to movement. It is possible that you also feel motivated to increase physical activity in everyday life. In addition, the activity sensor is provided with a sensor for measuring heart rate, which allows you to measure it in real time on your wrist. Thanks to the free "Medel Health Check" app, physical activity data can be easily transmitted from the activity sensor to the smartphone. In this way, you can view the recordings on your smartphone at any time, in a simple and intuitive manner.

In addition, you can receive notifications about messages (SMS, WhatsApp, Facebook Messenger) or incoming calls, and all this conveniently on your wrist. In this way,

everything is quickly under control and is always up to date.

MEDEL CARDIO WATCH with the "Medel Health Check" app offers the following functions:

- Date and time indication
- Display of the steps taken
- Display of daily distance travelled, calorie consumption, active time and progress of daily objective (percentage %)
- Representation of the individual phases of sleep
- Vibrating alarm clock
- Motion reminder
- Notification of calls and messages received on the smartphone (SMS, WhatsApp...)
- Search for smartphone
- Remote shooting with smartphone camera

System requirements:

- iOS $\geq$ Version 9.0 / Android™ $\geq$ Version 5.0

- Bluetooth® $\geq$ 4.0

The activity sensor uses Bluetooth® low energy technology, 2.402 - 2.480 GHz frequency band, maximum transmission power delivered in the <20 dBm frequency band, it is compatible with Bluetooth® 4.0 smartphones/tablets.

For the purpose of updating the technical data are subject to change without notice.

2. NOTES

WARNING

- Avoid contact with sunscreens or similar as they may damage plastic parts.
- Keep children away from packaging material. Danger of suffocation!
- Consult a doctor before starting a new training programme.
- This device is NOT a medical device. The measurement of the heart rate of the activity sensor is not suitable for diagnosing, treating, nursing or preventing diseases.
- The activity sensor contains magnets. Before using the activity sensor, consult your doctor if you are a pacemaker, defibrillator or other electronic system user
- The activity sensor is only provided with a splash-proof protection. It is however recommended NOT to use the activity sensor in the shower or while swimming.
- Protect the activity sensor from moisture.
- To protect the environment, the appliance must not be disposed of in household waste at the end of its use. Disposal must be carried out in the appropriate collection centres. Comply with local regulations in force regarding the disposal of materials. Dispose of the appliance according to the European Directive on Waste Electrical and Electronic Equipment (WEEE). For any clarifications, contact the municipal authorities responsible for disposal.

Warnings on the use of rechargeable batteries

- If the liquid from the rechargeable battery comes in contact with your skin and eyes, rinse the affected parts with water and consult your doctor.
- When a rechargeable battery is worn out, wear protective gloves and dispose of the device.
- Protect rechargeable batteries from excessive heat.
- **⚠ Risk of Explosion!** Do not dispose of rechargeable batteries in a fire.
- Do not break up, open or crush rechargeable batteries.
- Use only the chargers specified in the instructions for use.
- Rechargeable batteries need to be charged correctly before use. Follow the manufacturer's instructions and the information provided in these instructions for use to correctly charge the batteries.
- Before first commissioning, fully charge the rechargeable battery (see chapter 3.1)
- To achieve the longest possible life cycle, fully charge the rechargeable battery at least twice a year.

3. START-UP

To start-up the activity sensor, press and hold the button for 10 seconds. If the activity sensor does not turn on, the activity sensor's rechargeable battery is low. Charge the rechargeable battery as described in chapter 3.1: "Charging the battery".

1. Make sure your smartphone is connected to the Internet.
2. Activate the *Bluetooth*® on your smartphone.
3. Download the free "Medel Health Check" app from the Apple App Store (iOS) or from Google Play (Android™).
4. Start the "Medel Health Check" app and follow the instructions. In the app settings menu, select Medel CW and connect. After connection, you can set your personal data. Optionally, you can able/disable alarms, receive notifications for calls/messages or able/disable a motion reminder.

Tips to correctly connect the activity sensor to the smartphone:

- Position the smartphone as close as possible to the activity sensor.
- Keep a distance of 30 m from other *Bluetooth*® sensors/devices.

Depending on your smartphone settings, there may be additional costs for Internet connection or data transmission.

1. Charging the battery

Charge the activity sensor through the supplied USB cable, connecting it to your PC switched on or to a 5V dc 500 mA power supply. Be careful to position the cable in the correct position. During the charging process the battery symbol appears on the display.

To check the current battery charge status, press the button on the activity sensor. Generally a full battery charge requires about 3 hours. When the rechargeable battery is fully charged, the display shows the full battery symbol.

When fully charged, the rechargeable battery lasts approximately 7 days (with a heart rate measurement of one hour per day). The life of the rechargeable battery varies according to the lower or higher heartbeat measurement frequency.

Charge the rechargeable battery regularly to allow constant data recording.

 WARNING

DO NOT charge the activity sensor if the two charging contacts on the back of the device are wet. If necessary, gently dab the two charging contacts with a dry cloth. When a rechargeable battery is worn out, wear protective gloves and dispose of the device.

2. Wearing Medel CARDIO WATCH:

- Apply the activity sensor on the left wrist if you are right-handed or on the right wrist if you are left-handed.
- Apply the activity sensor on the wrist about 3-8 cm above the carpal bone.
- Make sure that the back of the light sensor is correctly in contact with the skin. The device must be tightened in such a way as to ensure that there is no light scattering of the sensor to ensure correct reading of heartbeats.

EN

If you often have cold or improperly blood supplied hands, it is recommended that you perform warm-up exercises before measuring your heart rate for a precise result.

3. Starting the sensor

Press the white O key to start Medel CARDIO WATCH:

If the activity sensor does not turn on, the activity sensor's rechargeable battery is low. Charge the rechargeable battery as described in chapter "Charging the battery".

4. FUNCTIONS

The activity sensor has the following functions: time, date and battery status, heart rate, pedometer, calorie consumption, distance travelled, active time, daily goal, alarm, motion reminder, sleep mode and messages/calls notification, smartphone search

EN

and remote shooting for the smartphone camera.
The values are reset every day at midnight. The count then starts again. To change function, press the activity sensor button.

4.1 Display of time / date / rechargeable battery


The activity sensor displays the current time, date and rechargeable battery status. When the activity sensor is synchronised with the smartphone, the current time and date of the smartphone are transmitted to the activity sensor.

4.2 Steps /movement

The activity sensor indicates movements in steps. Steps are recorded using a 3D sensor placed inside the activity sensor body.

4.3 Calorie consumption

The activity sensor calculates the calorie consumption based on the data specified in the app (height, age, gender, weight) and steps taken. When you wake up, the activity sensor already displays the number of calories consumed because the body consumes calories even in the rest condition. The appliance automatically considers the calories consumed in the rest condition.

4.4 Distance travelled

The activity sensor determines the length of the step through the height set in the app. The length of the step can also be customised. The activity sensor calculates

the distance travelled based on the length of the step and the number of steps. The displayed distance is only an approximate reference value and does not correspond exactly to the actual distance.

4.5 Active time 01:00

The active time indicates how long you have been on the move. The active time considers only the time during which one moves.

4.6 Progress of the daily goal 10%

You can configure a daily goal in the "Medel Health Check" app settings (e.g. 10000 steps). The progress of the set daily target (in %) can be viewed on the activity sensor display.

4.7 Alarm 19:45

You can set different alarm times in the "Medel Health Check" app. When the alarm time arrives, the activity sensor starts vibrating.

4.8 Messages / Notifications SMS John Smith WhatsApp John Smith

If push notifications are received on the smartphone (e.g. for SMS or WhatsApp), you can be warned on the activity sensor display. To do this, the activity sensor must be connected via Bluetooth® to the smartphone and the "Medel Health Check" app must be configured.

When the message arrives, the sender appears first on the display. If you press the


EN

button, part of the message is displayed. By pressing the button again you can continue reading the message. In case of longer messages, full contents are not displayed.

4.9 Calls  John Smith You can set the receipt of notifications on the activity sensor concerning calls in the "Medel Health Check" app settings. To do this, the activity sensor must be connected via *Bluetooth®* to the smartphone and the app must be configured.

4.10 Sleep mode To measure movements during sleep, activate the sleep mode of the device before going to sleep. To do this, press the button in the sleep mode indicator for 3 seconds. The display briefly shows 

The sleep mode is automatically activated after a second. When you wake up in the morning, press the button for approx. 3 seconds to reactivate the activity sensor from the sleep mode.

Optionally a time interval can be established in which the activity sensor automatically switches to the sleep mode and finally returns to the activity mode. You can also set a sleep reminder.

4.11 Activity alarm

You can set a motion reminder in the "Medel Health Check" app settings. A start time and an end time can be set for the interval in which the reminder is enabled, for example, from 9:00 am to 6:00 pm. In addition, an idle time can be set. If there is no movement during the set inactivity time, the activity sensor starts vibrating and appears on the display.  → .

4.12 Search for the smartphone

If it is impossible to find the smartphone, the activity sensor offers the possibility to make the smartphone sound and vibrate.

For that purpose:

- the two devices must be connected via *Bluetooth®*,
- the ringtones and the vibration must not be disabled on the smartphone,
- the distance between the two devices must be less than 5 metres,
- the "Medel Health Check" app must be opened on the smartphone.

1. Press the button on the activity sensor until the time or sleep mode appears on the display.

2. Press and hold the button on the activity sensor for 3 seconds. The display shows .
3. Now briefly press the button on the activity sensor. The display shows . The activity sensor now tries to find the smartphone. As soon as the smartphone is found, it will beep and start vibrating.
4. To stop the search for the smartphone and return to the main menu, press and hold the button on the activity sensor for 3 seconds. The display shows .

4.13 Remote shooting for the smartphone camera

With the activity sensor you can also trigger the smartphone camera. In this way, you can take pictures in a simple and practical way without having to hold your smartphone.

To enable the smartphone camera with the activity sensor:

- the two devices must be connected via *Bluetooth®*,
- the distance between the two devices must be less than 5 metres,
- the "Medel Health Check" app must be opened on the smartphone.

To trigger the smartphone camera through the activity sensor, proceed as follows:

1. Open the "Medel Health Check" app on your smartphone.
2. Press the button on the activity sensor until the time or sleep mode does NOT ap-

pear on the display.

3. Press and hold the button on the activity sensor for 3 seconds. The display shows

4. Press and hold the button on the activity sensor again for 3 seconds until the display shows

Now shortly press the button on the activity sensor. The display shows The smartphone's camera app starts.

6. Briefly press the button on the activity sensor again. The camera app on the smartphone only takes one picture. At each subsequent brief pressure of the button, the app takes a new picture.

4.14 Heart rate OFF

The activity sensor is provided with an optical sensor that measures the heart rate in real time. Heart rate measurement is only possible if the activity sensor is NOT in sleep mode.

1. To start measuring the heart rate, press the button until the display does not show

2. Press and hold the button for 3 seconds, the display shows and the activity sensor starts detecting the heart rate. As soon as the activity sensor has

detected a heart rate, the display shows  89 . If the activity sensor does not detect the heart rate, the display shows  OFF . In this case check if the light sensor on the back is correctly in contact with the skin and has no skin or dust particles. Then press the button to activate the heart rate measurement again.


3. To stop the heart rate measurement, press and hold the button for 5 seconds. The display shows  OFF

4.14 Display symbols

	Steps
	Heart rate
	Calorie count
	Covered distance (km)
	Duration of activity (h.-min.)
	Percentage Vs goal

	Target achievement
	Alarm function
	Camera control
	Find my phone
	Night/sleep mode

4.15 Overview of the functions

5. DATA TRANSMISSION TO THE SMARTPHONE

To transmit the recorded data from the activity sensor to the "Medel Health Check" app proceed as follows:

Enable the *Bluetooth®* on your smartphone.

2. Open the "Medel Health Check" app.

3. Connect the MEDEL CARDIO WATCH activity sensor to the app. Data is transmitted automatically to the smartphone.

6. CLEANING

Clean the activity sensor regularly with a damp cloth. Remove residues of skin and impurities from the strap and in particular from the plastic surfaces on the inside of the strap to avoid damage to materials. Also pay attention to the cleaning of the internal sensor for heartbeat detection, especially when using creams and/or cosmetics on the skin.

7. FREQUENTLY ASKED QUESTIONS

How is the activity sensor applied correctly?

Apply the activity sensor on the left wrist if you are right-handed or on the right wrist if you are left-handed. Thread one end of the strap through the opening of the other end.

What are the advantages of the activity sensor?

Movement is the fundamental condition for a healthy life. Physical activity helps keep

fit and reduces the risk of developing cardiovascular system diseases. Movement also helps to lose weight and increases well-being. Physical activity means any movement that involves energy consumption. Daily activities such as walking, cycling, climbing stairs and doing housework contribute greatly to movement. The activity sensor allows evaluating and objectively recording daily physical activity. It is possible that you also feel motivated to increase physical activity in everyday life. The ideal daily program should include 7,000 - 10,000 steps or 30 minutes of activity. (Source: according to WHO, Dr. Catrine Tudor-Locke, Dr. David R. Bassett jr.)

How does the activity sensor record the steps?

The activity sensor continuously records physical activity thanks to a built-in 3D sensor. All body movements are interpreted as steps. The indication of the route travelled may therefore differ from the actual journey.

How can the activity/sleep data be evaluated?

The *Bluetooth*® interface of the activity sensor allows transmitting stored data to a smartphone (AndroidTM and iOS) where they can be evaluated. Just download the free "Medel Health Check" app on your smartphone. The activity sensor uses *Bluetooth*® low energy technology, 2.4 GHz frequency band.

How are data synchronised with the app?

To transmit the recorded data from the activity sensor to the "Medel Health Check" app, follow the instructions in chapter "5. Data transmission to the smartphone".

EN

Data transmission to the smartphone does not work.

This can be due to several causes:

- Check if the smartphone/tablet is compatible with the "Medel Health Check" app.
- The activity sensor is in "Sleep mode". To exit the "Sleep mode", press and hold the button for 3 seconds.  Appears on the sensor display.
- The *Bluetooth*® has not been enabled on the smartphone. Enable the *Bluetooth*® in your smartphone settings.
- The activity sensor's rechargeable battery is low. If the activity sensor does not turn on, charge it with the supplied USB charging cable.
- Make sure you use the correct app. Data transmission is only possible with the "Medel Health Check" app.
- The activity sensor is connected to another smartphone. Interrupt the connection before connecting the activity sensor with another smartphone.
- **iPhone:** if you use an iPhone and the *Bluetooth*® connection is no longer established, interrupt the *Bluetooth*® connection first in the "Medel Health Check" app. Finally access the iPhone settings related to *Bluetooth*® and delete Medel CW from the list of *Bluetooth*® devices connected. Now connect Medel CW again in the "Medel Health Check" app via *Bluetooth*®.
- **Android™:** if you use an Android™ smartphone and the *Bluetooth*® connection is no longer established, remove the *Bluetooth*® connection in the "Medel Health Check" app and finally connect Medel CW again.

When should the activity sensor battery be charged?

Charge the activity sensor battery when, by pressing the button, on the display **19:45 New** appears. As soon as the display shows **19:45 New**, battery charging is limited to a maximum of 24 hours. Never fully discharge the battery pack. This is the only way to guarantee continuous data recording.

How must the activity sensor battery be charged correctly?

To charge the activity sensor correctly, follow the instructions in the "Charging the battery" chapter.

Can the activity sensor be worn even while riding a bicycle?

The activity sensor can also be worn while cycling, provided it is worn properly. For example, it is advisable to secure the sensor to the leg or in a narrow or side pocket of the trousers, do not wear it on the waist. However, the indication of the route travelled deviates from the actual route as the activity sensor interprets all body movements as steps.

How many days/how many nights can the activity sensor record?

The activity sensor can record 15 days and 15 nights. When all positions of the memory are occupied, the oldest recordings are automatically overwritten. It is therefore advisable to proceed with data transmission at least once a week.

How is the sleep mode enabled?

To do this, follow the instructions in the "Sleep mode" chapter.

How is the vibrating alarm clock enabled?

Set the desired alarm time in the "Medel Health Check" app. Synchronise the activity

EN

sensor with the app after entering/changing the alarm time.

Can the sensor be worn even in the shower or while swimming?

The activity sensor is protected against splashes of water, such as rain for example. In some cases, in the shower or during immersion of the activity sensor, it could happen that the prevailing pressure is too high and the activity sensor could be damaged. The activity sensor is not suitable for immersions.

How should the activity sensor be disposed?

To do this, follow the instructions in the "Disposal" chapter.

How must the activity sensor be cleaned?

From time to time carefully clean the activity sensor with a damp cloth.

What does the supply include?

The supply includes an activity sensor and a USB charging cable.

What do the symbols on the activity sensor and on the packaging mean?

To do this, see the chapter "Explanation of symbols".

8. DISPOSAL

Batteries must not be disposed of in household waste. They can contain harmful heavy metals and must be treated as special waste.

- The following symbols are shown on rechargeable batteries containing toxic substances:

Pb = leaded battery

Cd = battery containing cadmium

Hg = battery containing mercury

To protect the environment, the appliance must not be disposed of in household waste at the end of its use.

Disposal must be carried out in the appropriate collection centres. Dispose of the appliance according to the European Directive on Waste Electrical and Electronic Equipment (WEEE). For any clarifications, contact the municipal authorities responsible for disposal.

9. TECHNICAL SPECIFICATIONS

Weight: approx. 36 g (including cable)

- Suitable for wrist circumferences of approx. 14.5-19.5 cm
- Rechargeable battery: 150mAh/3.7V Polymer Lithium Ions
- Permitted operating conditions: 10-40° C
- Storage capacity: 15 days - Heart beats: 40-240 bpm
- Power supply: 5V DC, 500mA
- Water Proof: IPX4
- *Bluetooth*® low energy technology, 2.4GHz band frequency

We hereby guarantee that the product complies with European directive RED 2014/53/ EU. For further information contact the Medel International S.r.l. support service.

The *Bluetooth*® brand and the related logo are registered trademarks of *Bluetooth*® SIG, Inc.

The use of these trademarks by Medel International S.r.l. has been licensed. Other tra-

Subject to errors and changes

demarks and trade names are those of their respective owners.

10. EXPLANATION OF SYMBOLS

The following symbols are used in the Instructions for use and on the device:

	WARNING Reporting of risk of injuries or health hazards.		Manufacturer
	This product meets the requirements of current European and national directives.		Read the instructions for use
IPX4	Protected against solid objects over 1.0mm e.g. wires.	 5V DC, 500 mA	Power supply
	Disposal in accordance with the provisions of the EC Directive on waste electrical and electronic equipment (WEEE).		

11. WARRANTY

In the event of a breakdown, contact the Medel International S.r.l. Assistance centre. Never open the device. The appliance has no part inside it that can be repaired by unqualified personnel and does not require internal maintenance and/or lubrication.

Warranty conditions

The appliance is guaranteed for 3 years from the date of purchase against any original defect in materials or construction. The guarantee consists in the replacement of the defective parts at origin. The warranty does not cover the accessories supplied and the parts subject to normal wear and tear. The shipping costs of the appliance are charged to the user. Repairs outside the warranty conditions are charged to the user.

The warranty is void if the appliance has been tampered with, if the defect derives from improper use or if the damage is not attributable to the manufacturer (accidental fall, inaccurate transport, etc.). The guarantee does not entail any compensation for direct or indirect damages of any nature to persons or property during the period of malfunction of the product. The warranty is valid from the date of purchase of the product certified by the till receipt or the purchase invoice. In the event that the device does not resume proper operation, contact the MEDEL CUSTOMER SERVICE.

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