



PRODOTTI MEDICALI DAL 1996

La Sicurezza dagli Specialisti della Salute




D-Heart

D-Heart Portable ECG Device

Clinical Validations

Hospital Grade Multiple Leads Smartphone ECG



Medical grade ECG recording of 8-Leads ECG.

Accurate Morphologic evaluation of the ECG signal (DI, DII, DIII, aVR, aVL, aVR, V2, V5).

Innovative cable uncoiling mechanism designed to simplify and fasten the user experience.



Smartphone camera algorithm displays the electrodes position on the patient's chest with an accuracy of 0.5 cm*.

Personal medical folder with all medical data and exams history with the possibility to schedule exams.

24/7 telecardiology service that provides innovative triage report of the exam specifically designed for cardiac patients.

WHEN D-HEART CAN BE IMPORTANT FOR MY PATIENT



Arrhythmias assessment

Evaluation of palpitations

Rule out cardiac conditions as a cause of palpitations



Managing patients with cardiac conditions

Pre and Post Atrial Fibrillation procedure

Post Myocardial Infarction surveillance

Titration of newly introduced medications



Diagnosing Arrhythmias early in high risk patients

Post cardiac surgery

30-days after Transient Ischemic Attack



Regular Screening of patients with cardiac risk factors

Hypertension

Atypical Chest Pain

Cardiomyopathy




PRODOTTI MEDICALI DAL 1996

La Sicurezza dagli Specialisti della Salute



D-Heart

ECG Accuracy



Comparative analysis of multiple leads smartphone electrocardiograph (D-Heart®) versus standard 12-leads electrocardiograph in patients with Hypertrophic Cardiomyopathy

Niccolò Maurizi M.D. (1,2), **Carlo Fumagalli M.D. (1)**, Mattia Targetti M.D. (1), Silvia Passantino M.D. (1), Anna Arretini M.D. (1), Alessia Tomberli (1), Katia Baldini (1), Niccolò Marchionni M.D. (2), Franco Cecchi M.D. (3), Iacopo Olivotto M.D. (1)

(1) Unit Cardiomiopatie, Azienda Ospedaliero Universitaria di Careggi, Florence, Italy, (2) Department of Experimental and Clinical Medicine, University of Florence, Italy (3) Fondazione ARCARD, Florence, Italy.

Background. **M-Heath** technologies are revolutionizing cardiovascular (CV) medicine. However, a low-cost, user-friendly multi-lead smartphone electrocardiograph is still lacking. **D-Heart®** is a portable device that enables the acquisition of the ECG on multiple leads via 6 electrodes (3 peripheral, 3 augmented and two precordial leads -V2 and V5) which streams via Bluetooth to any smartphone. Both high and low-income settings may benefit from a low-cost/high-technology device.

Purpose. To assess the accuracy of D-Heart® recordings in the stratification of ECG morphological abnormalities, compared with 12-lead ECGs, in a cardiomyopathies outpatient clinic.

Methods. Consecutive patients (>18 years) with a diagnosis of Hypertrophic Cardiomyopathy (HCM, N=144, men=96) referred for outpatient control at a referral national institution for cardiomyopathies were enrolled from May to August 2017 (**Table 1**)

Table 1. Baseline Clinical characteristics	
	Overall
Demographics	
Population – N, (%)	144
Age	35±21
ICD/PM	16 (12%)
Echo Parameters	
LAD, mm	42±6
LVMWT, mm	21±5
Previous ACA	3 (1%)




Figure 1. D-Heart Smartphone ECG device

Results. Results of ECG abnormality and intervals measurement are summarized in **Figures 2, 3, and 4**. Agreement was obtained in 143/144 (99%) cases with D-Heart tracings and in 142/144 cases with 12-lead ECGs.





Figure 2. Agreement of ECG abnormalities (D-Heart vs 12-lead ECG) using a validated score




Figure 3. Distribution of R-E score points (D-Heart vs 12-lead ECG)

Figure 4. Examples of D-Heart ECGs recorded from HCM patients during the study



Figure 5. Comparison of PR and QRS intervals (Bland-Altman method, non-parametric approach) showed excellent concordance for D-Heart® measurements (95% limit of agreement -20 to +20 ms for PR and -10 to +10 ms for QRS).

Conclusions. **D-Heart®** proved effective and accurate, allowing stratification of ECG abnormalities comparable to the 12-lead electrocardiographs. These results open new perspectives for low-cost community cardiovascular screening programs in low-income settings or homecare delivery in high income countries.



PRODOTTI MEDICALI DAL 1996

La Sicurezza dagli Specialisti della Salute



Electrodes Placement Algorithm





PRODOTTI MEDICALI DAL 1996

La Sicurezza dagli Specialisti della Salute

Validation of a smartphone-camera based software for the identification of electrodes location on human chest

Niccolò Maurizi, MD^(1,2), Carlo Fumagalli, MD^(1,2), Mattia Targetti, MD⁽¹⁾, Silvia Passantino, MD⁽¹⁾, Anna Arretini, MD⁽¹⁾, Alessia Tomberli⁽¹⁾, Katia Baldini⁽¹⁾, Niccolò Marchionni, MD⁽²⁾, Iacopo Olivotto, MD⁽¹⁾, Franco Cecchi, MD⁽³⁾

⁽¹⁾ Cardiomyopathy Unit, Careggi University Hospital, Florence, Italy ⁽²⁾Department of Experimental and Clinical Medicine, University of Florence, Italy, ⁽³⁾ARCARD Foundation, Florence, Italy.

Background.

- Smartphone are producing profound changes in diagnostics.
- Clinically **actionable** data could now be generated locally by the patient.
- Multiple-lead **ECG** is a cost-effective, valuable and non-invasive test, but its use requires trained healthcare workers to correctly record the exam (due to specific position that electrodes).

Purpose.

To validate '**Mobile ECG on other**', a smartphone App AI algorithm that guides electrode self-placement via imaging processing technology by identifying electrode theoretical location on the patient's own chest, via the smartphone camera. The app is coupled with D-Heart®, a **multiple lead mobile ECG** designed for users with no medical background.

Methods. Methods are summarized in Figure 1 and 2.




Figure 1. Validation study Enrollment and analysis




Figure 2. A. Mobile on other AI electrode placement algorithm with D-Heart. B. D-Heart mobile ECG final placement

Results. Consecutive healthy volunteers (n=40, men=30) were enrolled from June to July 2017. Results are presented in Table 1, Table 2 and Figure 3.




Figure 3. Mobile on Other vs Nurse placement comparison

Table 1. Median displacement in the horizontal and vertical axis for each electrode in the modality D-Heart 'Mobile ECG on other'

	Δ Height	Δ Width
RA (cm)	0.3 [0.1-1] (0-3)	0.5 [0.3-1] (0-2)
LA (cm)	0.5 [0.2-1] (0-3)	0.3 [0.2-1] (0-3)
LL (cm)	1 [0.5-1] (0-2)	0.6 [0.5-1] (0-2)
V2 (cm)	0.5 [0.1-1] (0-3)	0.5 [0.1-0.9] (0-3)
V5 (cm)	0.5 [0.1-1] (0-3)	0.5 [0.2-1] (0-1.5)

Table 2. Number of electrodes misplaced by the algorithm in the modality D-Heart 'Mobile on other ECG'.

	Misplacement (N)	
RA	1	Six (15%) patients had >1 electrode misplaced.
LA	1	
LL	0	
V2	5	
V5	4	

Mean **time for ECG placement and recording** (defined as the time needed for the user to open the D-Heart soft case until the beginning of the ECG streaming from the D-Heart device to the smartphone) was **58 ± 12 seconds**.

Conclusions. 'Mobile ECG on other' software proved reliable and accurate for correct electrode self-placement, thereby opening new perspectives for accurate patient generated remote diagnostic tests.



PRODOTTI MEDICALI DAL 1996

La Sicurezza dagli Specialisti della Salute



D-Heart

Clinical Applications

ESPECIALLY VALUABLE FOR



Cardiologist



General
Practitioners



Nurse Practitioners



Remote Screening
Locations



Telemedicine
Diagnosis

Feasibility of cardiovascular screening in low-income settings using smartphone-based technologies

Carlo Fumagalli, M.D. (1)*, Niccolò Maurizi M.D. (1,2)*, Samuel O’Konu (3), Roberta Rughetti (3), Fulvio Avvantaggiato M.D. (4), Malamine Tamba (4), Mattia Targetti M.D. (1), Silvia Passantino M.D. (1), Anna Arretini M.D. (1), Alessia Tomberli (1), Katia Baldini (1), Fausto Barlocco PhD (1), Niccolò Marchionni M.D. (2), Franco Cecchi M.D. (5), Iacopo Olivotto M.D. (1)

(1) **Cardiomyopathy Unit, Careggi University Hospital, Florence, Italy**, (2) Department of Clinical and Experimental Medicine, University of Florence, Italy, (3) AMREF Health Africa, Nairobi, Kenya, (4) Comitato Pavia Asti Senegal, Ziguinchor, Senegal, (5) Fondazione ARCARD, Firenze, Italy. (*) Authors equally contributed to this project.

Introduction. Mobile health (**m-Health**) technologies revolutionizing cardiovascular (CV) medicine. However, more than **80% of heart disease related deaths** occur in **low income settings**, that are unable to take advantage of such innovations for healthcare awareness and prevention.

Purpose. To assess the feasibility of a pilot CV screening campaign with smartphone based m-health devices:
-D-Heart®, a validated low-cost 8-lead electrocardiogram (199€)
-iHealth® blood pressure (BP) recorder (~70€)

Methods. A total of **231 patients** were enrolled in a two-day screening at 4 rural dispensaries, 2 in Kitui District (**Kenya**) and 2 in Ziguinchor District (**Senegal**).







Figure 1. Blood Pressure measurement results **Figure 2.** ECG Screening results



Results. Clinical characteristics of the 231 patients are summarized in **Table 1**, **Figure 1** and **Figure 2**.

Table 1. Baseline characteristics of the screened population

	Overall	Men	Women	p
Demographics				
Population – N, (%)	231	121 (52)	110 (48)	
Age	36±21	41±17	32±12	<0.01
BMI	24±2.9	22.3±1.5	27.2±1.1	<0.01
CV risk Factors				
Smoking Hx – N, (%)	78 (34)	60 (50)	18 (16)	<0.01
Alcohol (>2.5 U/day) – N, (%)	92 (40%)	67 (55)	37 (33)	<0.01
HTN Hx – N, (%)	11(5)	5 (4)	6 (5)	ns
Diabetes Mellitus	15	7 (6)	8 (7)	ns
Previous AMI	3 (1)	2 (2)	1 (1)	ns

71% never had BP measured before

91% never had ECG recording before

Figure 3. Samples of D-Heart ECGs recorded

Cost-Effectiveness Analysis. Taking into account:

- device price
- consumables
- salaries
- USB solar powered technology
- visit time (6±2 minutes).

€1.10/patient

- **€0.80 for community health worker**
- **€0.30 for consumables**

Conclusions. D-Heart® ECG screening combined with smartphone BP measurement proved **efficient** and **cost-effective**. This should encourage to develop **low-cost/high-technology** community-based **CV screening programmes** in low-income settings.

Disclosures. Dr. Fumagalli was awarded a research grant from D-Heart Srl in December 2017

