



## Statement

Our Coronavirus Ag Rapid Test Cassettes (Swab) and Influenza & COVID-19 Ag Combo Rapid Test Cassette (Swab) are capable of detect the mutated COVID-19 501.v2 south African Variant . Our test cassettes mainly detect the nucleoprotein of the COVID-19 501.v2 south African Variant. The third amino acid residue of N gene of COVID-19 was not located in antigenic determinant region, and the antigenicity of N gene was not affected after mutation. Therefore,

Our Coronavirus Ag Rapid Test Cassettes and Influenza & COVID-19 Ag Combo Rapid Test Cassette (Swab) can be used to detect the new variant of COVID-19 denominated ( COVID-19 501.v2 south African Variant )which is currently Mutated in the South African . However, our antigen test cannot be used for the differential diagnosis of COVID-19 and mutated COVID-19.

Name: Victoria Zhang  
Position: Sales Manager

